

Guide to spheroid formation using Verified for Spheroids™ Human Hepatocytes

Use this protocol for spheroid formation with our Verified for Spheroids™ Human Hepatocytes to support your next steps towards more physiologically relevant liver cell culture systems.

Introduction

The liver is the body's main site of metabolism of most small molecule drugs and consequently, *in vitro* models that accurately replicate liver function are highly desirable. It has recently been shown that primary hepatocytes, the main functional cells of the liver, can self-assemble into small spheroids when cells are placed into low-adherence plates or hanging liquid drops.¹ The formation of these spheroids improves the *in-vivo*-like response and increases the life span of the hepatocytes. In this Technical Note, we describe the development of a simple spheroid formation protocol using our Verified for Spheroids™ Human Hepatocytes, to support researchers desire to implement hepatocyte spheroids as a routine model for *in vitro* metabolism, disease modeling, and toxicity testing.

Materials

1. Primary Verified for Spheroids™ Human Hepatocytes (Lonza, cat. no. HUCPI or HUCPG)
2. HCM™ Bullet™ Kit (Lonza, cat. no. CC-3198)
3. Fetal Bovine Serum (FBS, Hyclone, cat. no. SH30071.03)
4. 1M HEPES Buffer (Lonza, cat. no. 12509079)
5. Human Cryopreserved Hepatocyte Thawing Medium (Lonza, cat. no. MCHT50)
6. 96-well Corning U-bottom ultra-low attachment (ULA, Sigma, cat. no. CLS7007)

Protocol

1. To formulate HCM™ Hepatocyte Culture Medium, transfer the contents of the HCM™ SingleQuots™ Kit to HBM™ Basal Medium with a pipette, and rinse each vial with medium. Store at 4°C for up to 1 month.
2. Make complete Spheroid Formation Medium (SFM) by adding the following supplements to HCM™ Medium: FBS (final 20%) and HEPES Buffer (final 25mM). Store at 4°C for up to 1 month.

Note: We tested the effect of various levels of serum addition to the media on spheroid formation. Results indicate that spheroids require a minimum of 10% serum to form, but 20% serum produces the most optimal spheroids (Figure 1).

3. Thaw Verified for Spheroids™ Human Hepatocytes in thawing medium and centrifuge at speeds and times described in our Protocol. <https://lonza.picturepark.com/Go/ACI6hFB1/D/29887/1>

Note:

Each batch of Lonza's HUCPI and HUCPG are characterized for spheroid formation. Look for the Verified for Spheroids™ Specification in batch certificates or visit product web pages for up-to-date inventory.




Figure 1

At least 10% serum is required for spheroid formation. 3,000 primary human hepatocytes seeded in ULA U-bottom plates with serum at concentrations ranging from 0 – 20%. Optimal results are obtained when 20% serum is included in spheroid formation medium.

4. Resuspend cells in 3 – 8 mL SFM and count cells using trypan blue and a hemocytometer.
5. Adjust cells to a final concentration of 1.0×10^4 – 3.0×10^4 cells/mL in SFM depending on the size of spheroid desired.

Note: We tested seeding between 1,000 – 5,000 cells per well. While all seeding densities produced spheroids, increasing the number of cells/well resulted in larger spheroids with cores beyond 500 μ M. This enlarged size increases the risk for necrotic cores and therefore, we recommend using less than 3,000 cells/well when forming spheroids (Figure 2).




Figure 2

Spheroids of different size. 1,000 to 5,000 primary human hepatocytes cultured in ULA U-bottom plates on day 7 after second media change.

Spheroids reach size of 500 μ M when 3,000 cells are used for spheroid formation.

6. To begin spheroid formation, add 100 μ L cell suspension to wells of a U-bottom ULA 96-well plate and place in incubator at 37°C 5% CO₂. Incubate undisturbed for 120 – 168 hours.

Note: For extended cultures of spheroids, we recommend filling all outer edge wells with medium only to assist with evaporation issues. Any disturbance of the plates or mishandling during incubation can negatively affect spheroid formation.

7. Once spheroids have formed, remove 50% of medium and add fresh complete HCM™ Medium every 2 to 3 days. More frequent media changes can facilitate faster removal of remaining serum used during spheroid formation.

Our Verified for Spheroids™ Human Hepatocytes are plateable human hepatocytes that are already qualified as general purpose (cat. no. HUCPG) or Interaction Qualified (cat. no. HUCPI) and further tested for forming a tight spheroid by 7 days in culture. In our extended studies, it was observed that there is some donor-to-donor variability in both the formation rate and longevity of hepatocyte spheroids.

Additional studies can be found in a separate White Paper, <https://lonza.picturepark.com/Go/dfhpfzvf/D/34835/1> that show certain hepatocyte lots can last in spheroid culture up to 28 days with continued high-level Cytochrome P450 activity. These efforts support the potential use of our Verified for Spheroids™ Human Hepatocytes for extended toxicology and DMPK studies and reduce the need for researchers to test different donors for spheroid formation.

References

¹ Bell CC, Hendriks DF, Moro SM, Ellis E, Walsh J, Renblom A, Puigvert LF, Dankers AC, Jacobs F, Snoeys J, Sison-Young RL. Characterization of primary human hepatocyte spheroids as a model system for drug-induced liver injury, liver function and disease. Scientific reports. 2016 May 4;6:25187.

Contact information

North America

Customer Service: +1 800 638 8174 (toll free)

Fax: +1 301 845 8338

order.us@lonza.com

Scientific Support: +1 800 521 0390 (toll free)

scientific.support@lonza.com

Europe

Customer Service: +32 87 321 611

order.europe@lonza.com

Scientific Support: +32 87 321 611

scientific.support.eu@lonza.com

International

Contact your local Lonza distributor

Customer Service: +1 301 898 7025

scientific.support@lonza.com

Lonza Walkersville, Inc. – Walkersville, MD 21793

For research use only. Not for use in diagnostic procedures.

All trademarks belong to Lonza or its affiliates or to their respective third party owners. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights.

©2019 Lonza. All rights reserved.

CD-DS035 01/19

www.lonza.com/pdh