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## FAQ

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Standard Vertical Systems	PAGEr™ ゲル
PAGEr™ Minigel Chamber	9 cm × 10 cm 10 cm × 10 cm*
Bio-Rad® MiniPROTEAN® II, MiniPROTEAN® 3, Mini-PROTEAN® Tetra, Mini-PROTEAN® Dodeca™ and Ready Gel® Cell Systems Reverse the inner core gasket so the flat side faces outward.	9 cm × 10 cm
Novex® XCell SureLock® Mini-Cell or XCell II Request the spacer for the XCell SureLock®, Mini-Cell Chamber from Scientific Support, [Cat. No. 59900].	9 cm × 10 cm 10 cm × 10 cm*
FisherBiotech® Vertical Minigel FBVE121, Owl Separations Systems Wolverine™ P82 Chamber comes with 2 sets of wedges. Use the thinner wedges for the PAGEr™ Gold ゲル.	10 cm × 10 cm
FisherBiotech® Vertical Minigel FB-VE101, Owl Separations Systems Penguin™ Model P8DS Request adaptor for these chambers from Scientific Support, [Cat. No. 59902].	10 cm × 10 cm
Hofer® Mighty Small™ (SE250) Replace the buffer chamber with a 'Deep lower buffer chamber for the SE260', order number 80-6148-78, from GE Healthcare.	9 cm × 10 cm 10 cm × 10 cm*
Daiichi 2, ISS chambers To run one gel: Place one 10 × 10 cm cassette on wedge side of chamber. Use suitable buffer dam on the other side. Use regular Daiichi/ISS wedges. To run two gels: Widen the hole on the yellow port of the inner core. Replace the long arm wedges with modified wedges. This chamber modification and new wedges are available from Scientific Support.	10 cm × 10 cm
Hofer® Mighty Small™ (SE260)	9 cm × 10 cm 10 cm × 10 cm
EC 120 Mini Vertical Gel System	9 cm × 10 cm 10 cm × 10 cm
Biometra® Mini V Chamber	9 cm × 10 cm
CBS Scientific MGV System, (10 cm × 8 cm units)	9 cm × 10 cm
Sigma-Aldrich Mini Techware (11.3 cm × 10 cm units)	10 cm × 10 cm
Zaxis System 2000	10 cm × 10 cm
Hofer® Mini VE	10 cm × 10 cm

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\*Omit for native proteins.

For best results use Lonza AccuGENE™ Electrophoresis Buffers

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NOTE: Limits are based on optimal detection methods for each stain.

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Nitrocellulose	PVDF	Nylon
Hydrophobic binding	Hydrophobic binding	Hydrophobic and electrostatic binding
General purpose membrane	SDS tolerant	Stable if baked
Low background	High background	High background
Low strength	High strength	High strength
Becomes brittle if baked	Suitable for protein sequencing	Least suitable for Western transfer

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