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
図1. The general path of stable integration.

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Culture System	Advantage	Application
Batch culture — polyclonal	Fast, useful for cells which do not grow in single cell culture	Overexpression, protein expression systems [e.g., for basic research]
Limiting dilution — monoclonal	Defined cell clones	Studies of gene function, protein production [e.g., for therapeutic applications]

In a batch culture system, a mixed population of drug resistant cells is selected on plates or in flasks and can be used directly for experimental analysis. During a limiting dilution procedure, cells are usually diluted and selected e.g., in a 96-well plate for outgrowth of cell clones or single colony growth. Subsequently, colonies can be picked and used to generate monoclonal cell lines.

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






図2. Higher integration rates in difficult-to-transfect cell lines using Nucleofection™. Jurkat cells were transfected using either Nucleofection™ (2 µg DNA) or lipofection Reagent L (0.7 µg DNA) according to the respective manufacturer's instructions. 24 hours after transfection, cells were plated on a 96-well plate containing culture medium supplemented with G418 for selection of stably-transfected cells. 30 days after plating, cells were analyzed for clonal outgrowth [Integration frequency = number of resistant clones per number of living cells seeded]. Due to toxic effects, lipofection with 10 µg circular DNA could not be performed.

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	Procedure Outline	Important Information
Step 1 	Design experiment and choose cell type, expression vector and transfection method.	Make sure that transfection method and expression vector are suitable for your cell type.
Step 2 	Determine appropriate cell number per well (only for limiting dilution) and G418 concentration.	Cells differ in their susceptibility to G418. The active concentration of stock G418 can vary from batch to batch.
Step 3 	Transfect expression vector into cells.	Amount of expression vector per experiment is dependent on transfection method and cell type.
Step 4 	Plate transfected cells and cultivate cells in medium without G418.	Do not add G418 to culture medium immediately after transfection as this may drastically increase mortality.
Step 5 	Dilute cells into culture plates and start selection 24 – 48 hours post-transfection. Feed every 2 – 3 days (for batch culture) or 10 days (for limiting dilution) with selection medium.	Choose culture conditions (batch culture, limiting dilution) depending on your experimental design. Refreshed selection medium is important to avoid false positive cells.
Step 6 	Analyze stably transfected cells.	Make sure that the chosen assay is suitable for your application.

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- Grimm S. [2004]. The art and design of genetic screens: mammalian culture cells. Nature Rev Gen, 5: 179-189.
- Hacein-Bey-Abina, S., Von Kalle, C., Schmidt, M., McCormack, M. P., Wulffraat, N., Leboulch, P., Lim, A., Osborne, C. S., Pawliuk, R., Morillon, E., Sorensen, R., Forster, A., Fraser, P., Cohen, J. I., de Saint Basile, G., Alexander, I., Wintergerst, U., Frebourg, T., Aurias, A., Stoppa-Lyonnet, D., Romana, S., Radford-Weiss, I., Gross, F., Valensi, F., Delabesse, E., Macintyre, E., Sigaux, F., Soulier, J., Leiva, L. E., Wissler, M., Prinz, C., Rabbitts, T. H., Le Deist, F., Fischer, A., Cavazzana-Calvo, M. [2003]. LM02-associated clonal T cell proliferation in two patients after gene therapy for SCID-X1. Science, 302: 415-419.
- Southern, P. J., Berg, P. [1982] Transformation of mammalian cells to antibiotic resistance with a bacterial gene under control of the SV40 early region promoter.
- J Mol Appl Genet, 1: 327-341.
- Wurm, F. M. [2004]. Production of recombinant protein therapeutics in cultivated mammalian cells. Nature Biotechnol, 22: 1393-1398.