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
Figure 2. Chondrogenic Differentiation hMSC were differentiated with the hMSC Differentiation BulletKit™ – chondrogenic (Lonza, PT-3003) in the presence or absence (undifferentiated control) of 10ng/ml TGF-β3 (Lonza, PT-4124). Chondrogenic differentiation was monitored on the indicated days using Alcian blue 8GX staining [Sigma, A3157]. Absorption was measured with the FLUOstar Omega (BMG Labtech) at 610 nm [data is expressed as mean ± SD for 3 replicates]. All lots underwent chondrogenic differentiation as measured by glycosaminoglycan production compared to undifferentiated controls. Lot-to-lot dependent variability was observed.


Figure 3. Osteogenic Differentiation hMSC were differentiated with the hMSC differentiation BulletKit™ – osteogenic (Lonza, PT-3002) according to manufacturer's instructions with exception of the seeding density [1,85E+03 cells/cm²]. Osteogenic differentiation was monitored on the indicated days with the Osteolmage™ Mineralization Assay (Lonza, PA-1503). Fluorescence intensity was measured with the Fluoroskan Ascent [Labsystems; excitation 485 nm/emission: 520 nm; 293 data points/well; data is expressed as mean ± SD for 1 – 3 replicates]. All donors underwent osteogenic differentiation as measured by hydroxyapatite deposition of cells. Lot-to-lot dependent variability was observed with respect to total amounts of deposited hydroxyapatite per well as well as the time required for detectable hydroxyapatite deposition.



Figure 4. Adipogenic Differentiation hMSC were differentiated with the hMSC Differentiation BulletKit™ – adipogenic (Lonza, PT-3004) according to manufacturer's instructions. All donors underwent adipogenic differentiation as measured by the accumulation of intracellular lipid droplets within the cells. The percentage of differentiated cells and/or the total lipid content per cell varied between different lots.

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