

Dulbecco's Phosphate
Buffered Saline
w/o Calcium
or Magnesium

17-512

Lonza

Description	CAS #	Chemical Formula	Concentration		Molarity	
			g/L	mg/L	M	mM
Potassium Chloride	7447-40-7	KCl	0.200	200.000	2.683E-03	2.683
Potassium Phosphate Monobasic Anhydrous	7778-77-0	KH ₂ PO ₄	0.200	200.000	1.470E-03	1.470
Sodium Chloride	7647-14-5	NaCl	8.000	8.000E+03	0.137	136.893
Sodium Phosphate Dibasic-7-Hydrate	7782-85-6	Na ₂ HPO ₄ • 7H ₂ O	2.160	2.160E+03	8.058E-03	8.058