

Dulbecco's Phosphate
Buffered Saline
w/o Calcium
or Magnesium (10X)

Lonza

17-515

Description	CAS #	Chemical Formula	Concentration		Molarity	
			g/L	mg/L	M	mM
Potassium Chloride	7447-40-7	KCl	2.000	2000.000	2.683E-02	26.828
Potassium Phosphate Monobasic Anhydrous	7778-77-0	KH ₂ PO ₄	2.000	2000.000	1.470E-02	14.695
Sodium Chloride	7647-14-5	NaCl	80.000	80000.000	1.369	1368.925
Sodium Phosphate Dibasic-7-Hydrate	7782-85-6	Na ₂ HPO ₄ • 7H ₂ O	21.600	2.160E+04	8.058E-02	80.576