

Phosphate Buffered Saline
w/o Calcium, Magnesium,
or Phenol Red

17-516



Description	CAS #	Chemical Formula	Concentration		Molarity	
			g/L	mg/L	M	mM
Sodium Chloride	7647-14-5	NaCl	9.000	9.000E+03	0.154	154.004
Sodium Phosphate, Dibasic Anhydrous	7558-79-4	Na ₂ HPO ₄	0.795	795.000	5.599E-03	5.599
Potassium Phosphate Monobasic Anhydrous	7778-77-0	KH ₂ PO ₄	0.144	144.000	1.058E-03	1.058