

Phosphate Buffered Saline
w/o Calcium, Magnesium, 17-517
or Phenol Red (10X)



Description	CAS #	Chemical Formula	Concentration		Molarity	
			g/L	mg/L	M	mM
Potassium Phosphate Monobasic Anhydrous	7778-77-0	KH ₂ PO ₄	1.440	1.440E+03	1.058E-02	10.580
Sodium Chloride	7647-14-5	NaCl	90.000	9.000E+04	1.540	1.540E+03
Sodium Phosphate, Dibasic Anhydrous	7558-79-4	Na ₂ HPO ₄	7.950	7.950E+03	5.599E-02	55.986