

Phosphate
Buffered Saline-EDTA

BE02-017



Description	CAS #	Chemical Formula	Concentration		Molarity	
			g/L	mg/L	mM	uM
Na ₂ HPO ₄ anhydrous	7558-79-4	Na ₂ HPO ₄	0.713	713.246	5.023	5.023E+03
KH ₂ PO ₄	7778-77-0	KH ₂ PO ₄	0.232	231.971	1.704	1.704E+03
NaCl	7647-14-5	NaCl	8.766	8.766E+03	150.000	1.500E+05
KCl	7447-40-7	KCl	0.280	280.000	3.756	3.756E+03
EDTA Na ₂ ·2H ₂ O	6381-92-6	C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ · 2H ₂ O	0.460	460.000	0.878	878.214
Lithium Chloride	7447-41-8	LiCl	0.427	426.700	10.066	1.007E+04