

Versene (EDTA) 0.02%

17-711

Lonza

Description	CAS #	Chemical Formula	Concentration		Molarity	
			g/L	mg/L	mM	uM
Potassium Chloride	7447-40-7	KCl	0.200	200.000	2.683	2.683E+03
Potassium Phosphate Monobasic Anhydrous	7778-77-0	KH ₂ PO ₄	0.200	200.000	1.470	1.470E+03
Sodium Chloride	7647-14-5	NaCl	8.000	8.000E+03	136.893	1.369E+05
Sodium Phosphate Dibasic-7-Hydrate	7782-85-6	Na ₂ HPO ₄ • 7H ₂ O	2.170	2.170E+03	8.095	8.095E+03
EDTA (Ethylenediaminetetraacetic Acid Disodium Dihydrate)	6381-92-6	C ₁₀ H ₁₄ N ₂ O ₈ Na ₂ •2H ₂ O	0.200	200.000	0.537	537.288