

## Grace's Insect Medium

---

### Product Use

Grace's Insect Medium is used to cultivate insect cell lines and organs for the propagation of entomopathic viruses.

### Specifications

	<b>Sterility</b>	<b>pH</b>	<b>Osmolality (mOsm)</b>
04-457	Negative	5.9-6.2	340-380
04-501	Negative	FIO	FIO
04-649	Negative	FIO	FIO

  

	<b>Cell Growth Generation</b>	<b>Endotoxin (EU/ml)</b>
04-457	Pass	FIO
04-501	≥75% of control	FIO
04-649	≥75% of control	FIO

### Storage

2°C to 8°C

### Product Use Statement

**THESE PRODUCTS ARE FOR RESEARCH USE ONLY.** Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* procedures.

### Ordering Information

Catalog Number	Description	Size
04-457F	Grace's Insect Medium without lactalbumin hydrolysate, yeastolate, gentamicin or FBS	500 ml
04-501Q	Grace's Insect Medium complete with 10% heat inactivated FBS, yeastolate, lactalbumin hydrolysate and gentamicin	1 L
04-649F	Grace's Insect Medium with lactalbumin hydrolysate, yeastolate and gentamicin without FBS	500 ml